

HYTO GmbH

Hochfeld 242

5721 Piesendorf

Austria

Bischofshofen,
14.12.2017

Appraisal 17318691

Processability of Hyto Toilet System

Following inspection of the Hyto toilet system product for processability based on ISO 17664 (see Inspection Report 17318611 of 30.11.2017), we have evaluated the product as follows:

The Hyto toilet system can be safely cleaned through treatment in a cleaning and disinfection unit with a standard cleaning program and a conventional cleaning agent for medical products (e.g. mediclean forte). A cleaning performance > 99.9% was verified in testing.

The Hyto toilet system can be safely disinfected in a cleaning and disinfection unit with a standard thermal disinfection program for "ward utensils" (target value 80°C/1 minute and A0 value of 60 s).

The Hyto toilet system can be safely sterilised through treatment in a steam sterilizer at a temperature of 134°C.


Dr. Arno Sorger
Technical Director, Inspector

To: Client